



GRS Genomic DNA Kit – Blood & Cultured Cells –
#GK02.0100
(FOR RESEARCH ONLY)

SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM BACTERIA

GRAM-NEGATIVE BACTERIA

- 1) Transfer up to 1×10^9 cultured bacterial cells to a 1,5-ml microcentrifuge tube, centrifuge for 1 minute at 14.000g-16.000g and discard the supernatant
- 2) Resuspend the bacterial pellet in **200 μ l of Buffer BC2** and incubate at room temperature for 5 minutes. Proceed with the lysis step of the cultured cell protocol (page 6).

GRAM-POSITIVE BACTERIA

- 1) Transfer up to 1×10^9 cultured bacterial cells to a 1,5-ml microcentrifuge tube, centrifuge for 1 minute at 14.000g-16.000g and discard the supernatant
- 2) Resuspend the bacterial pellet in **200 μ l Lysozyme Buffer¹** and incubate at room temperature for 10 minutes. During incubation, invert the tube regularly. Proceed with the lysis step of the cultured cell protocol (page 6).

¹ Lysozyme Buffer = 20mg/ml Lysozyme; 20mM Tris-HCl pH 8.0; 2mM EDTA; 1% Triton X-100
(prepare fresh immediately prior to use)