

Xpert cDNA Synthesis Mastermix

GK81.0100 (100 rxn)
(FOR RESEARCH ONLY)

Product Description

Xpert Reverse Transcriptase (RNase H-) is a RNA-dependent DNA polymerase suitable for cDNA synthesis from long RNA templates. Xpert RTase has been optimized to perform under high temperatures (45°C-55°C), which facilitates the removal of secondary mRNA structures associated with high GC content. Together with the lack of RNase H activity, which ensures minimization of template degradation during long incubation times, Xpert RTase enables the preparation of long full-length cDNAs (up to 15kb).

The Xpert cDNA Mastermix provides an efficient and fast method for the synthesis of high quality first strand cDNA from purified poly(A)+ mRNA or total RNA templates. The Xpert cDNA Mastermix (2X) is an optimized reaction mix containing a balanced concentration of oligo(dT) and random hexamer primers, dNTPs, and RNase inhibitor. Note that, since the mastermix already contains primers, it cannot be used with gene-specific primers. First strand cDNA can be directly used as template in PCR.

Kit Contents

Component	Volume (100 rxns)
Xpert Reverse Transcriptase (200U/μl)	100 μl
Mastermix (2X)	1,2 ml
RNase-free water	2 ml

Quality Control

Functionally tested by cDNA synthesis and PCR amplification in comparison with previous batches. Absence of endonuclease-, exonuclease-, and nicking activity is verified on a lot-to-lot basis.

Storage

Store all components at -20°C for up to 1 year.

Protocol

1. Mix the following components in a RNase-free microtube:

Component	Volume
template RNA	1 ng - 2 µg total RNA or 1 pg – 2 ng poly(A)+ RNA
Mastermix (2X)	10 µl
RNase free water	up to 19 µl

2. Using a thermocycler or thermoblock, heat for 5 minutes at 65°C (in order to remove possible secondary RNA structures), then place on ice for 2 minutes.
3. Add the following components to the mixture:

Component	Volume
Xpert RTase (200U/µl)	1,0 µl

4. Mix thoroughly and then centrifuge briefly
5. Incubate at 25°C for 10min
6. Using a thermocycler or thermoblock, heat the microtube for
 - a. 15 min at 50°C (for usage in qPCR) or
 - b. 50 min at 50°C (for usage in PCR)
7. Inactivate Enzyme by heating for 5 min at 85°C
8. Either use cDNA immediately as template in qPCR/PCR or store at -20°C.

Notes

1. The RNA sample should be completely free of contaminating genomic DNA.
2. There is no need to purify poly(A)+ RNA from total RNA when using oligo(dT)₂₀ primers, however, doing so, may improve yield and overall purity of the final product.
3. One can prepare a no RT control by taking a small aliquot of the Xpert RTase and inactivate the enzyme by incubating at 85°C for 5-10 minutes, prior to adding in step 3.

GRiSP Research Solutions

Rua Alfredo Allen, 455
4200-135 Porto
Portugal
www.grisp.pt | info@grisp.pt