

Trypsin/EDTA (1x solution)

#GTC02.0100
(for research only)

| | |
|------------------------|---|
| Product: | Filter sterilized solution of 0.05% Trypsin from porcine pancreas in Dulbecco's PBS (pH 7.0-7.6) with EDTA, without Ca ²⁺ , without Mg ²⁺ , and without phenol red. Ready-to-use cell detachment solution for breaking down cell adhesion structures. |
| Quantity: | 100ml |
| Appearance: | Clear colourless solution. |
| Storage: | -20°C for up to 2 years. |
| Specifications: | Activity: > 500 BAEE U/ml pH: 7.0–7.6 Osmolality: 250-320 mOsm/kg Sterility: sterile Parvovirus test : negative Performance: Cell Detachment |

Typical trypsinization protocol

(depends on cell line and should be adapted conform)

Prior to start, pre-warm all solutions (Trypsin/EDTA (1x) and PBS (without Ca²⁺ and without Mg²⁺)) to 37°C by placing in a water bath. The entire protocol should be carried out in a laminar flow hood, using proper aseptic techniques.

1. Carefully aspirate off all the culture media from the flask, without disturbing the cell layer or letting the cell layer dry. Proceed immediately with washing the cells with PBS (rinse with 5 ml for a T25 flask and with 10 ml for a T75 flask).
 2. Immediately after, add sufficient Trypsin/EDTA (1x) solution to cover the cells. For a T25 flask, 0,5ml* should be enough. Incubate the flask at 37°C for 2-3 minutes in the cell culture incubator. Check cell morphology visually (microscope) to verify if cell have rounded, if not continue incubation. The solution in the flask will appear cloudy. Incubation time should be kept at a minimum and overexposure should be avoided, as Trypsin can cause cellular damage. The required time depends among others on cell type, culture age, cell density, and the serum concentration in the growth medium).
 3. Once cells are rounded and detached, beat the flask against the palm of your hand to loosen any remaining attached cells
 4. Wash out all the cells from the surface by pipetting the fresh complete cell culture medium (5ml) all over the surface. Gently disperse the cells to break cell clumps. Take a sample to determine the viable cell density, and add aliquots of detached cells to fresh culture media in new flasks.
- *) When using more, e.g. 2-4ml, then in step 4, you must add more serum containing medium to inhibit trypsin after digestion has been completed, or neutralize excess of trypsin by adding trypsin inhibitor.

GRiSP Research Solutions

Rua Alfredo Allen, 455
4200-135 Porto
Portugal

www.grisp.pt | info@grisp.pt