



**GRS Genomic DNA Kit – Tissue –**  
**#GK03.0100**  
(FOR RESEARCH ONLY)

---

**SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM AMNIOTIC FLUID**

- 1) Transfer approximately 10ml (up to 15ml) of amniotic fluid to a 15 ml centrifuge tube (not provided) and harvest the cells by centrifugation for 3 minutes at 14.000g-16.000g. Discard the supernatant.
- 2) Resuspend the pellet in 200µl of **Buffer BC2** and transfer to a 1,5-ml microcentrifuge tube. Add 10µl of **Proteinase K** (10mg/ml), mix by shaking vigorously, and incubate at **60°C** for 30 minutes. During incubation, invert the tube regularly.
- 3) Add 200µl of **Buffer TC1**, mix by shaking vigorously, and incubate at **60°C** for at least 20 minutes. During incubation, invert the tube regularly. [At this time, preheat the Elution Buffer in a 60°C water bath].
- 4) [*Optional step; If RNA-free DNA is required*] Allow the mixture to cool to room temperature and add 5µl of **RNase A (10mg/ml)**, mix by shaking vigorously and incubate for 5 minutes at room temperature.
- 5) Add 200µl of absolute ethanol to the lysate and mix immediately by shaking vigorously for 10 seconds. In case precipitate appears, break it up by pipetting.
- 6) Place the **genomic DNA mini spin column** in a 2-ml collection tube and transfer the sample mixture (including any precipitate if present) to the column.
- 7) Centrifuge at 14.000g-16.000g for 1 minute. Discard the collection tube containing the flow-through and place the genomic DNA mini spin column in a new collection tube.
- 8) Add 400µl of **Wash Buffer 1** and centrifuge at 14.000g-16.000g for 30 seconds. Discard the flow-through and place the genomic DNA mini spin column back in the collection tube. Add 600µl of **Wash Buffer 2\*** and centrifuge at 14.000g-16.000g for 30 seconds (\*check if ethanol is added first time the kit is used; see Notes on page 2).
- 9) Discard the flow-through and place the genomic DNA mini spin column back in the collection tube and centrifuge for another 3 minutes at 14.000g-16.000g to dry the matrix of the column.
- 10) Transfer the spin column to a new 1,5-ml microcentrifuge tube and pipet 100µl **preheated Elution Buffer** directly to the centre of the spin column without touching the membrane. Incubate at **37°C** (incubator) for 10 minutes.  
**Notes:** Instead of Elution Buffer, DNA can also be eluted with TE or water; pH should be 8.0-8.5. Standard elution volume is 100µl. To increase concentration elute with 30-50µl. To increase yield, elute with 200µl.
- 11) Centrifuge for 1 minute at 14.000g-16.000g to elute purified genomic DNA. Discard the spin column and use DNA immediately or store at -20°C