



GRS Genomic DNA Kit – Tissue –
#GK03.0100
(FOR RESEARCH ONLY)

SUGGESTED PROTOCOL FOR DNA PURIFICATION FROM BUCCAL SWAB

- 1) Add 500µl of **Buffer BC2** and 20µl of **Proteinase K** (10mg/ml) to a 1,5-ml microcentrifuge tube. Place the buccal swab into the tube and incubate at **60°C** for 10 minutes.
- 2) Discard the swab and add 500µl of **Buffer TC1** to the lysate. Mix immediately by shaking vigorously. Incubate at **60°C** for 10 minutes. [At this time preheat the elution buffer to 60°C for the DNA elution step].
- 3) [*Optional step; If RNA-free DNA is required*] Allow the mixture to cool to room temperature and add 4µl of **RNase A (10mg/ml)**, mix by shaking vigorously and incubate for 5 minutes at room temperature.
- 4) Add 500µl of absolute ethanol to the lysate and immediately mix by shaking vigorously.
- 5) Place the **genomic DNA mini spin column** in a 2-ml collection tube and transfer 700µl of the sample mixture. Centrifuge at 14.000g-16.000g for 1 minute.
- 6) Discard the flow-through and repeat step 5 with the remaining lysate mixture.
- 7) Proceed with the wash step of the Tissue Protocol on page 3 (step 9).